

Development of high density (600K) chicken genotyping array

Almas Gheyas
Roslin Institute, University of Edinburgh
almas.gheyas@roslin.ed.ac.uk

15th Jan, 2012 Plant and Animal Genomes (PAG) XX Conference















Steps in array development



Major steps:

- Re-sequencing of many chickens and detection of SNPs
- Selection of SNPs for validation
- Validation of selected SNPs
- Final selection of 600K SNPs for array



Fig: Affymetrix Axiom™ genotyping array
















Re-sequencing



- Illumina high throughput sequencing
- 10-15 samples pooled within line
- Average depth of coverage: 7-17x per line















SNP detection



- Sequence reads were aligned to chicken reference genome (Gallus_gallus-4.0)
- SNP detection performed using Samtools (v .0.1.7a)
- SNP detection done by
 - Within line analyses: 78M SNPs
 - All lines together: 139M SNPs
















SNP selection




















Validation of 1.8M SNPs



- 1.8M SNPs placed in three Affymetrix arrays for validation
- 288 individuals genotyped:
 - 3 broiler lines
 - 9 layer lines
 - 26 outgroup individuals

















SNP validation criteria



Good genotype cluster resolution

(included in the validated SNP set)

- SNP must be polymorphic with ≥3 copies of minor allele
- SNPs must be robust
 - Genotype call rate ≥98%
 - Good cluster separation
 - Reproducibility
- SNPs have stable Mendelian inheritance
- SNPs must meet performance criteria across all lines

Polymorphic and robust SNPs	1,187,482	65%
Stable Mendelian inheritance	1,176,808	99%















Selection of final 600k



- Backbone was created by taking all validated SNPs within exonic gene regions (n=21,534)
- SNPs with extreme breach of HWE (P<0.00001) were removed.
- LD between adjacent pairs of SNPs were calculated (Golden Helix SVS) and the average broiler and layer LDs were used to decide the optimum balance of broiler and layer SNPs (3:2).
- Broiler and layer SNPs were selected to get as even spacing as possible in terms of map distance.
















Annotation of 600K SNPs



















600K SNP distribution






Fig 3: Distribution of SNPs across chromosomes in unit physical (mb) and map distance (cM)

















Inter-marker spacing







Mean gap size is $1,748 \pm 5,274$ bases
















Principal component analysis

PROSLIN

PCA analysis was performed to see if the markers can group the closely related individuals

















Uses of high density SNP array ROSLIN

- Genomic selection
- Selection signature analysis
- Genome wide association studies
- For fine mapping of QTLs















Acknowledgements



Funding: BBSRC/Defra LINK grant; Aviagen Ltd, Affymetrix Ltd; Hyline International; German Federal Ministry of Education and Research

Roslin Institute: <u>David Burt</u>; <u>John A. Woolliams</u>; <u>Chris Haley</u>; Almas Gheyas; Clarissa Boschiero; Frances Turner; Andy Law; Le Yu; Peter Kaiser; Paul Hocking

Aviagen: Kellie A. Watson; Andreas Kranis

Hyline: Janet E. Fulton

ARK genomics: Richard Talbot; Frances Turner; Sarah Smith; Alison Downing

Affymerix: Fiona Brew; Ali Pirani; Lucy Raynold; J. Montgomery; J. Bart

Synbreed: Henner Simianer; Ruedi Fries; Rudolf Preisinger; T.M. Strom; Steffen Weigend; Klaus Meyer; George Haberer













